

# BRUNO

# Scooter Lift Jr.®



◆ The Bruno Scooter-Lift Jr. Model ASL-400 Vehicle Lift ends the struggle of manually trying to lift and store a scooter in the trunk or hatch opening of your vehicle. The Scooter-Lift Jr. is designed to help raise and store three-wheel scooters into most auto trunks, station wagons or hatchbacks, as long as the opening is large enough to store the partially disassembled scooter and compact lift. It gives you the convenience you need to transport your scooter and the freedom to travel when and where you want.

◆ The Scooter-Lift Jr. is electrically powered from the vehicle's battery. Simply press a button on the controller and the Scooter-Lift Jr. raises/lowers your scooter, effortlessly. It can also be adapted for operation in both station wagons and hatchbacks. This versatility can provide a significant benefit when a new vehicle purchase is not in your budget.

◆ The Scooter-Lift Jr. is also available with the added option of a convenient spring-assist that further reduces the effort needed to raise and lower the lift arm before/after use.

## OPTIONAL BRUNO POW'RPAC



- ◆ Easy to use!
- ◆ Convenient!
- ◆ Works in most minivans, SUVs, station wagons & sedans!

The complete, independent power supply for your automotive lift!

Pow'rPac easily installs in your vehicle – just plug into the vehicle's power outlet!

◆ **EXCLUSIVE!** Bruno's R<sup>3</sup> Lift Belt is 50% stronger! Count on Bruno's R<sup>3</sup> Lift Belt for ease of use, longer life and safer operation – *another Industry First!*

◆ Expert installation and complete, after-sale service is provided by Bruno's worldwide network of dealers.

◆ Only Bruno has over 25 different interior/exterior scooter, wheelchair and power chair lifts for nearly any sedan, minivan, sport utility vehicle (SUV), full-size van, pickup truck and the new cross-over vehicles.

◆ All Bruno Lifts are eligible for reimbursement under the terms of Mobility Programs offered by GM, Saturn, DaimlerChrysler/Dodge, Ford, Toyota, Lexus, Volvo and Volkswagen.



# Using the *Scooter Lift Jr.™* is easy!

1



Once alongside the auto's rear bumper, little effort is needed to raise the lift arm that latches automatically when placed in the **UP** position.

2



You can swing the arm out up to 90°, from your auto, using little effort. Then, hook the *lift claw* to the exclusive Bruno docking device.

3



The scooter is ready to be raised. Just push the **UP** button and raise your scooter until it clears the 'trunk lip'.

4



The scooter must be positioned over the trunk opening. The controls to lower the scooter are conveniently located on the lift arm.

5



Releasing the lift arm for storage requires that the latch be raised and the lift arm be lowered to the trunk floor.

6



When the operation is completed, your scooter fits snugly in your auto trunk. The seat can be stowed in the trunk or inside the car.

**Bruno lifts are eligible for reimbursement under the terms of Mobility Programs offered by GM, Saturn, DaimlerChrysler/Dodge, Ford, Toyota, Lexus and Volkswagen.**



The Lift weight rating is the **MAXIMUM** load value. Factors inherent to specific applications including, but not limited to, vehicle structure and condition, operating environment, and customer ability may reduce the effective weight rating.

**Represented By:**



Bruno Independent Living Aids, Inc., 1780 Executive Drive,  
P.O. Box 84, Oconomowoc, WI 53066

1-262-567-4990 ♦ FAX: 1-262-953-5501

[www.bruno.com](http://www.bruno.com)